

General Specifications

Nominal Wavelength	755nm
Pulse Duration	750ps
Pulse Energy	200mJ
Electrical Power	200-240 V~,4.5 kVA, 50/60Hz, single phase
Electrical Service Requirement	30Amp

Handpiece Specifications

Zoom Handpiece <i>(tattoo removal)</i>	
Spot Size	2.0-6.0mm
Zoom Handpiece <i>(tattoo removal)</i>	
Spot Size	6.0, 8.8 10.0mm



For the power of PicoSure contact your Cynosure consultant today.

Call 800-886-2966. Visit www.picosure.com

REFERENCES:

1. Harris Interactive, 2012
2. Kilmer, et.al. <http://www.emedicine.com/derm/topic563.htm>. Accessed March 2009.
3. Patient's Guide survey, published May 22, 2012. <http://blog.patientsguide.com/2012/05/study-tattoo-removal-climbs-32-in-competitive-job-market/>
4. ABC News, May 2012
5. http://esciencenews.com/articles/2008/07/21/study_examines_motivations_tattoo_removal
6. Successful and Rapid Treatment of Blue and Green Tattoo Pigment With a Novel Picosecond Laser, Jeremy A. Brauer, MD; Kavitha K. Reddy, MD; Robert Anolik, MD; Elliot T. Weiss, MD; Julie K. Karen, MD; Elizabeth K. Hale, MD; Lori A. Brightman, MD; Leonard Bernstein, MD; Roy G. Geronemus, MD



Removing tattoos
just got faster

PicoSure

